

*a*RUNA

HYALURONIC ACID MTS KIT

ARUNA



HYALURONIC ACID MTS KIT

USING PREMIUM HOME CARE PRODUCTS, 3,600 MICRONEEDLES PENETRATE THE SKIN AND DELIVER COLLAGEN, SOOTHING, AND MOISTURIZING INGREDIENTS DEEP INTO THE SKIN, INCREASING SKIN ELASTICITY, REDUCING WRINKLES, AND IMPROVING MICROCIRCULATION.

HYALURONIC ACID MTS KIT

Let's do aesthetic care and home care at the same time!

#Once of Rolling Only

#Moisture Hydration Effect

#Effective Ingredient Absorption



01

Dual functionality for whitening & wrinkle improvement

Contains niacinamide and adenosine to help with skin whitening and wrinkle improvement.

02

Moisturizing effect

Contains ingredients that help moisturize the skin, keeping the skin hydrated.

03

Soothing effect on the skin

Contains madecassoside to help soothe the skin.

04

Contains peptides

Contains peptides to help with skin elasticity and conditioning.

05

Contains 12 types of amino acids

Contains 12 types of amino acids to provide nourishment and support skin health.

aRUNA

HYALURONIC ACID MTS KIT

Recommended for:

I recommend it to people like this

Premium home care!

Delivers active ingredients deep into the skin for radiant and elastic skin!

#Firm Skin

#Fast Absorption

#Healthy Skin Care



Those who find it difficult to visit a skincare shop and feel burdened by maintenance costs

Those who feel the limitations of the effects of topical skincare products

Those who concern about gradually losing skin elasticity

Those who want their skincare products to be more effectively absorbed

Those who want to replenish moisture in rough skin

aRUNA

HYALURONIC ACID MTS KIT

Product introduction



ARUNA MTS ROLLER

Aruna MTS Roller

Derma rollers stimulate the skin, maximizing the absorption of skincare products.

#Enhanced Effects on Skin

#Clean Skin

#Micro Circulation Care

#Nutrient Absorption



ARUN HA AMPOULE

Aruna HA Ampoule

Contains PDRN to aid in skin regeneration and help cultivate healthy skin.

#PDRN ampoule

#Pink Vitamin Included

#Firm Skin

#12 Types Of Amino Acids



CORE CELL AMPOULE

Aruna Core Cell Ampoule

Contains over 250 types of growth factors to moisturize the skin and promote flexibility.

#Hydrated Skin

#Moisture Soothing

#Elasticity Care

#Smooth Skin Texture

aRUNA

HYALURONIC ACID MTS KIT

Home care program

Premium Home Care' completed once a week

5-week Home Care Program




MTS Roller X 1ea + MTS Refill X 5ea
+
HA Ampoule X 5EA + Core Cell Ampoule X 5EA

※ THE ROLLER IS RECOMMENDED FOR USE ONCE A WEEK.
AFTER USING THE ROLLER ONCE, USING A MASK PACK OR
SERUM DURING THE 6-DAY REST PERIOD HELPS WITH SOOTHING.

The start of solving Skin concerns is **Absorption!**

WITH 3,600 MICRO-TIPS,
DEEPLY INTO THE SKIN!



WHEN USING HANDS
COTTON PADS



WHEN USING
THE MTS ROLLER



The needles create pathways in the skin,
helping the product penetrate deep into
the skin and between wrinkles.

aRUNA

HYALURONIC ACID MTS KIT

Micro Tip

Micro Tips for effective absorption
Of active ingredients into the skin.




The Micro Tips, made with active ingredients
such as Sodium Hyaluronate,
2 types of peptides, and Madecassoside,
dissolve to deliver these active ingredients deeply
and quickly into the skin.

aRUNA

HYALURONIC ACID MTS KIT

MTS roller point ingredient



DUAL FUNCTIONALITY FOR WHITENING & WRINKLE IMPROVEMENT

Niacinamide
Adenosine




MOISTURIZING EFFECT

Sodium Hyaluronate,
Tremella Fuciformis Polysaccharide



SKIN SOOTHING

Madecassoside



2 TYPES OF PEPTIDE

Acetyl Hexapeptide-8,
Acetyl Octapeptide-3

※ Subject to raw material characteristics



ARUNA AMPOULE

HA AMPOULE & CORE CELL AMPOULE

Ampoule Point



**DUAL FUNCTIONALITY
FOR WHITENING &
WRINKLE IMPROVEMENT**

Niacinamide
Adenosine




MOISTURIZING EFFECT

Sodium Hyaluronate,
Tremella Fuciformis
Polysaccharide



SKIN SOOTHING

Madecassoside



2 TYPES OF PEPTIDE

Acetyl Hexapeptide-8,
Acetyl Octapeptide-3

※ Subject to raw material characteristics



ARUNA HA AMPOULE

POINT INGREDIENTS

CONTAINS SALMON DNA, PDRN

PDRN, similar to human DNA, helps with skin care, moisturization, and nutrition, contributing to firm, healthy skin.

CONTAINS PINK VITAMIN (B12)

Contains cyanocobalamin to help with skin elasticity and moisturization, and to quickly soothe problematic skin.

SKIN EXFOLIATION EFFECT

Contains citric acid to regulate pH, effectively removing dead skin cells and helping to manage a healthy stratum corneum.

※ Subject to raw material characteristics



ARUNA CORE CELL AMPOULE

POINT INGREDIENTS

Contains Human Umbilical Mesenchymal Stem CellConditioned Media Extract

As a culture supernatant, which is secreted during stem cell culture, is rich in over 250 growth factors and cytokines. And the growth factors in the skin help regulate skin balance, aiding in skin elasticity management.



Micro Tips for effective absorption of active ingredients into the skin.



Check out our domestic and international certifications for the MTS roller, ensuring safe and reliable use!

※ Subject to raw material characteristics

ARUNA

HYALURONIC ACID MTS KIT

How to use

Aruna ampoule Usage instructions

With Aruna ampoule, deeply deliver nutrition

* Open the outer transparent container of the ampoule,
then open the lid and attach the dropper.
Use it effectively with the roller to apply the desired amount.



- ① Designed in an appropriate size to fit securely in the silicone puff's groove to prevent it from falling out.
- ② Made with large grooves for easy gripping.
- ③ The handle length is extended for a more stable grip.




Sterilization treatment

Clean and safely
sterilized



Melting chip protection

Protects the contents
with special foil



Melting chip support

Enhances silicone puff
thickness and elasticity



aRUNA

HYALURONIC ACID MTS KIT

MTS kit how to use

TO DEEPLY DELIVER NUTRITION INTO THE SKIN WITH AMPOULE & ROLLER

STEP.1



Open one portion of the roller and attach it to the top of the dedicated handle.

STEP.2



After cleansing, massage with the roller on skin without moisture or oil for 5-10 minutes.

STEP.3



Apply the ampoule to the area massaged with the roller.

STEP.4



Repeat the roller massage once more for effective absorption of the serum's active ingredients

ARUNA

HYALURONIC ACID MTS KIT

Q&A



Precautions for use

Please do not use the microchip roller for any other purpose than its original intended use.

It is advisable to refrain from using it on areas with wounds or cuts.

If your skin is excessively dry, using the roller may cause more irritation.

The needles melt upon contact with moisture, so ensure to dry the skin with a towel after cleansing, before using the needles.

A desiccant is included to maintain humidity please do not use it for any other purposes.

Q. How often should I use it?

This contains 5 units, so using it once a week allows for one month of home care. We recommend using it 1-2 times a week and allowing a minimum of 3 days for skin recovery.

Q. Can I use the roller alone?

We recommend using it with an ampoule or serum to effectively deliver the active ingredients contained in the roller.

Q. When using it on the skin, isn't it stimulating?

There may be slight stimulation during use, but Please rest assured that Aruna MTS Solution Kit uses a roller made entirely of active ingredients, providing a sensation of deep absorption into the skin.