



IRIS

BLUE TONING™

Beyond
Your Youth™

Bluetoning gives
your youth and more...



THE **WORLD BEST** TECHNOLOGIES SHOW:

1. Easy Installation
2. Easy Maintenance
3. Easy Monitoring
4. Easy Warranty
5. Easy Education

WHAT YOU GET? YOU WILL GET:



PERFECT TONING TREATMENTS

HIGH PEAK POWER

ENOUGH POWER CONTROL

TRUE 15 Hz

MAX ENERGY 1600 / 2000mJ

VARIOUS PTP MODES

EASY INSTALLATION

EASY MAINTENANCE

BEFORE SERVICE

REMOTE ASSISTANTS


SYSTEM MONITORING

SMARTPHONE ALERT

TRUE MAXIMUM ENERGY

The maximum energy test with a 9 ϕ rod and super hive flat lens at the end of the handpiece unlike other machines

1600mJ / 2000mJ



Advance Laser Optic Parts




TRUE 15 Hz

Power control technology supports true 15Hz all the time as much as possible without any energy variation.




True 15 Hz

PTP BLUETONING TECHNOLOGY (6:4)



Bluetoning controls the 1st and 2nd Q energy rate for each skin types and user preferences in the world to give dramatic treatment results.



THE WORLD FIRST CENTRAL MONITORING SYSTEMS

1. Auto alert system
2. Central monitoring system
3. Remote service system
4. Auto upgrading system



EASY HANDLING

Soft buttons and physical button with a jog dial provide user conveniences with various short cut menus



Jog Dial



Short Cut Keys

INTUITIVE GRAPHIC INTERFACE (IGI)

Users find out how to use intuitively without manual according to Iris' graphic user interface (IGI)



Basic GUI



Fractional GUI

EASY MAINTENANCE & WARRANTY

Iris has a slot type circuit electronic power system like a normal PC and it has auto calibration mode for handpieces




Auto-Calibration



Circuit Slots

NO NOISE WITHOUT WATER TANK & FAN

Iris doesn't have a plastic water tank with the patented technology for cooling system. It doesn't have loud noise.



No Water Tank



No Fan

OPTIONAL HANDPIECES



KTP 532 Fractional



585nm Gold



1064 Fractional



650nm Ruby



1064 7 Fixed Collimated

INDICATIONS

Iris has proven technology and power treatments with toning, whitening, peeling and also wide range of other indications for all skin types

- Melasma
- Freckle
- Inflammatory Acne
- PTP Bluetoning
- Abnom
- Soft Peel (&SLP Mode)
- Tattoo Removal
- Nevus of Ota
- Cafe-au-lait

SPECIFICATIONS

Device Name	Iris			
Laser	Q-switched Nd:YAG			
Beam Profile	Super Hive Flat mode (top hat mode)			
Mode	1064	532	SLP	PTP Blue Toning
Pulsed Energy[mJ]	1600 / 2000	500	2200~3000	1600 / 2000
Wavelength[nm]	1064	532	1064	1064
Pulse Width	5~10 ns	5~10 ns	300 μ s	5~10 ns
Spot[mm]/ Repetition[Hz]	2~10 / 1~10	2~10 / 1~10	2~10 / 1~10	2~7 / 10~15 8~9 / 10~12 10 / 10
Display	10.2" wide panel and jog dial smart panel			
Auto Management	Ethernet, WIFI, WCDMA			
Optional Part	7mm single collimated / 7x7 fractional handpeice (1064nm, 532nm)			
Memory	User programmable 8 memories to save and load parameters			
Noise Level	50dB (patent cooling system)			
Spot Size	Auto sensing zoom collimated handpeice (2 ~ 10mm)			
Lamp Life Time	10,000,000 shot or more			
Electrical Requirements	220~230VAC, 15A , 50Hz or 60Hz			
Dimensions(mm)/Weight	295(W) x 740(L) x 973(H) , 70Kg			



IRIS

PTP BLUETONING THERAPY



BLUECORE COMPANY CO., LTD.

#1203 Centum jungang-ro 48, Haeundae-gu, Busan 48059, Republic of KOREA

Tel. +82-51-747-4318 | Fax. +82-51-747-4319

Bluecorecompany.com | info@bluecorecompany.com

©2017 All Rights Reserved. Bluecore Company Co., Ltd. Blue Toning® is a registered trademark of Bluecore Company Co., Ltd. IRIS™ is a trademark of Bluecore Company Co., Ltd. Specifications are subject to change without notice. **IRIS is FDA Cleared.**